

SegVator[®] LeanSeat[™]





SegVator[®] is a registered trademark and LeanSeat[™] is a trademark of SegVator[®], LLC. Segway[®] is a registered trademark of Segway[®], Inc. SegVator[®] is not manufactured or endorsed by Segway[®], Inc.





The Segway[®] Personal Transporter (PT) is part of the mobility equation – the SegVator[®] LeanSeat[™] is another part







- High and Low LeanSeat[™] versions available
- Interchangeable between
 Gen 1 and Gen 2
 Segway[®] PT's Each
 LeanSeat[™] version will
 work with both
- Surprisingly comfortable





- Split "walk through"
 flip-down design
 allows you to use your
 Segway® PT
 conventionally with
 either LeanSeat™
 version
- Easy to quickly install and remove





- LeanSeat[™] utilizes an internal rod and "V" plate system to activate the platform sensors while seated
- Rugged, light-weight
 & compact with simple
 bolt-on design





LeanSeat[™] users of
 either High or Low
 version must be able to
 quickly stand
 unassisted and
 dismount the Segway[®]
 PT in the unlikely
 event of a mechanical
 or safety shutdown







- High LeanSeat[™] version useful for anyone with standing duration issues to extend glide time
 - For use when stationary to rest and "take a leg break" without having to dismount the Segway® PT
 - Not intended for extended or aggressive gliding while seated









- Low LeanSeat[™] version can be used for extended gliding
 - Must be able to fully sit and return to a standing position unassisted
 - Removable gel cushion available for added comfort





SegVator® proceeds
 benefit Disability
 Rights Advocates for
 Technology (DRAFT)
 & the Segs4Vets
 program –
 details available at
 www.Segs4Vets.org







Contact your authorized Segway[®] dealer for further information regarding the SegVator[®] LeanSeat[™]



Expanding Mobility with Universal Design

SegVator, LLC

11200 W. Silver Spring Rd.

Milwaukee, WI USA 53225-3198

414/461-9100

414/615-2197 Fax

www.segvator.com